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
Technical Data	105 II	205 II	305 II	315 II	325 II	425 II
wavelength (nm)	139 x 103 = 14283 532					
cw output power (mW)	> 10	> 20	> 50	> 100	> 150	> 400
transversal mode/ roundness of beam	TEM ₀₀ (>95%) / < 1.1 : 1.0					
TEM ₀₀ beam diameter, typical (mm) (1/e ²) at beam waist	0.32 *					
TEM ₀₀ beam divergence (mrad) (1/e ²)	2.2 *					
beam pointing stability (% of beam divergence)	< 1.0					
stability of output power over 8 hrs (%)	< ±2 (typical 0.5)					
noise (<10Hz to >1GHz)(% rms)	< 0.5 (typical 0.1)					
polarization	linear 100:1, vertical					
inherent linewidth	< 10 kHz					
jitter	<± 500 kHz					
drift	< 100 MHz/h (ΔT <±1 °C)					
operating voltage	110/220V AC ±10%					
power consumption	< 80W (typ.40W)			<200W (typ.100W)		<300W (typ.150W)
ambient temperature range for operation	0-35°C (32-95°F)**					
laser head heat sink temperature for conductive cooling	0-55°C (32-130°F)					
dimensions of laser head mm (inches), weight	275x55x51 (10.8x2.6x2.0) 1.2 kg			335x89x80 (13.2x3.5x3.1) 6.5 kg		
dimensions of power supply mm (inches), weight	215x105x170 (8.5x4.1x6.7) 4 kg			440x260x130 (17.3x10.2x5.1) 13 kg		

Vibration, shock (IEC 68)

Design and technical data are subject to change without notice.


* factory provided collimating optics on request

** with factory provided or other adequate heat sink

optional: micro-processor control unit with RS 232 interface (215 x 105 x 45 mm)

warm up time 5 min

Replacing tubes with semiconductors


ADLAS GmbH & Co.KG
Seelandstraße 67
23569 Lübeck, Germany
Tel. ...49-451-3909300
Fax....49-451-3909399



ADLAS, Inc.
636 Great Road
Stow, MA 01775, USA
Tel.: (508) 897-0800
Fax: (508) 897-0811

100/200/300-series

laser head


heat sink



320/420-series

laser head


heat sink

532 nm cw, Nd:YAG, single longitudinal mode


Diode Pumped Solid State Lasers

ADLAS manufactures diode pumped solid state lasers for a wide range of scientific and OEM applications. In power performance and beam quality our lasers demonstrate considerable advantages over other lasers. The principle choice of a laser diode pump over conventional lamp pumping is based on solid reasons: the narrow emission bandwidth of the diode laser matches the absorption band of the Nd:YAG crystal (see diagram to the right). Temperature stabilization of the carefully selected diode laser guarantees that the ideal wavelength of 808 nm is constantly maintained.


absorption band of Nd:YAG laser

ADLAS lasers feature a special pumping configuration (see schematic below) which offers several advantages: low pump power loss through direct coupling, high efficiency through mode matching of the pump and Nd:YAG laser beams, and a compact, rugged design. The high power, low beam quality pump light of the diode laser is efficiently converted (typ. 35%) to the diffraction limited output of the solid state laser. Additional noteworthy features are the low power requirement and the long, maintenance-free lifetime.

emission spectrum of diode laser

Diode pumped solid state lasers from ADLAS are cost effective in procurement and operation. Their stability and compactness make them perfectly suited to many purposes, and, as a result of their near-perfect beam quality, destined for maximum performance applications.

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
2.2(2) 2/28

ADLAS: Effective Laser Solutions

Diode pumped solid state lasers manufactured by ADLAS have provided solutions for many applications:

- spectroscopy
- metrology
- holography
- interferometry
- wafer inspection/-processing
- micro-material processing
- range finding
- medicine (diagnostic/therapeutic)
- printing
- light scattering (static/dynamic)
- particle counting
- photoluminescence
- signal transmission

and many more. Commercially viable solutions for further applications are under development.

Biotechnology

Biotechnology, a highly dynamic scientific-technological sector, requires new tools of extraordinary quality and precision. Typical applications for ADLAS lasers are "optical tweezers" (trap, hold, and move micron size particles), confocal microscopy, fluorescence excitation, cell counting and cell sorting.

Good reasons to choose ADLAS lasers:

- excellent focussability due to a near-perfect Gaussian laser beam
- high beam pointing stability
- very stable output power level
- easily adapted as add-on (e.g. to microscopes)
- compact design
- operable in any position


TEM₀₀ beam of an ADLAS laser

General Specifications

Beam quality	TEM ₀₀ ($M^2 < 1.05$)
polarization	linear > 100 : 1
noise (10 Hz - 1 GHz)	< 0.5 % rms (typ. < 0.1 % rms)
output power stability (8 h time period)	< ±2 % (typ. < ±0.5 %)
repetition rate of q-switched laser	from single shot to 50 kHz* (infrared QS lasers can be switched to cw operation)
beam pointing stability	< 1 % of beam divergence

*(reduced pulse energy for repetition rates >1 kHz)
Detailed specification sheets are available.

General features of ADLAS lasers

- compact, rugged design
- no high voltage
- sealed cavity design
- high efficiency
- air cooling (no liquids)
- extraordinary beam quality
- highly stable output power
- long coherence length
- TEM₀₀ mode
- high signal-to-noise ratio
- low power requirement
- maintenance-free
- minimal operating costs
- long lifetime

Test and Measurement


Test and measurement technology applications call for more specifically tailored laser products than ever before. ADLAS currently offers about 40 different approved and tested standard lasers, for example the q-switched lasers with ns pulses for range finding and the frequency doubled cw lasers with extremely long coherence length for non-destructive testing (holography, ESPI and shearography). Further applications are in spectroscopy and interferometry.


Good reasons to choose ADLAS lasers:

- very long coherence length
- short warm-up time
- excellent temperature stability
- rugged construction
- quiet operation
- air cooling
- compact design
- adjustable output power

Courtesy of Fraunhofer-Institute IPT, Aachen*

Product Overview

cw Lasers IR

1064 nm	50 mW - 3000mW
1313 nm	40 mW - 800 mW

cw Lasers VIS, SLM

532 nm	10 mW - 400 mW
--------	----------------

cw Lasers IR, SLM

1064 nm	20 mW - 600 mW
---------	----------------

Q-Switched Lasers

	models up to:	
1313 nm	>100 μ J	<70ns
1064 nm	>200 μ J	<20ns
1047 nm	>350 μ J	<20ns
532 nm	>60 μ J	<17ns
523 nm	>120 μ J	<17ns
349 nm	>40 μ J	<13ns
262 nm	>20 μ J	<13ns

Options


- C-mount
- microprocessor control with RS 232 interface (532 nm cw models only)
- adjustable output power
- polarization preserving SM fiber pigtail
- 30 dB isolator integrated in laser head (1313 nm cw models only)

Communication

Applications in the communications industry, cable television (CATV), broadband transmission, antenna remoting etc. provide a growing market for solid state lasers. ADLAS meets this demand by supplying lasers with integral optical isolators and permanently coupled fiberoptic pigtails. A cost-effective transmission system utilizes an ADLAS Nd:YLF solid state laser (1313nm) in conjunction with an external modulator. This type of system configuration, with its high output power, substitutes the need for several directly modulated systems at a lower cost.

Good reasons to choose ADLAS lasers:

- long lifetime (design for >10 years available)
- single line operation
- distortion free signal
- TEM₀₀ mode
- high signal-to-noise ratio
- polarization preserving coupling


Printing

In the printing industry, computer typesetting and raster image processing call for constantly shorter production times and ever increasing resolutions. The latest generation of printing plates are now able to be imaged directly. The printing environment offers many challenges to the components involved (vibrations, large temperature fluctuations). Our factory sealed resonator and its robust construction uniquely qualify ADLAS lasers to meet just such special requirements.

Good reasons to choose ADLAS lasers:

- low noise for clearly better exposure results
- no periodic fluctuations, smooth transitions
- excellent beam pointing stability
- simple cooling requirements
- low power consumption/heat load
- high output power for shorter exposure cycles


Material Processing

ADLAS q-switched lasers are well suited for material processing applications where micron precision is needed. Resistor trimming, flat panel display repair, silicon wafer marking, ASIC chip production and tuning of quartz crystals are just a few examples of ADLAS lasers on the production line. Coincidentally, every second quartz watch manufactured in Europe has been tuned with an ADLAS laser. As most of these processes take place in a Clean Room environment the benefits of maintenance-free, no-water operation, of the compact design and of an extremely small footprint become obvious.

Good reasons to choose ADLAS lasers:

- superior beam quality
- short, high energy pulses
- excellent pulse-to-pulse stability
- ease of operation

ADLAS
advanced design lasers

High Quality Products

QUALITY is of utmost importance to us at ADLAS. Quality management is an integral part of all functions within the company. Early failure recognition and avoidance procedures are being followed from the conceptual phase all the way to the finished product.

Design

Product design balances stringent component quality standards with cost effectiveness. Specific customer requirements can be addressed during this phase.

Development

Our engineers constantly evaluate components, products and applications to assure the quality standards demanded by the user. This process also leads to ongoing improvement of existing products.

Documentation

To assure the product integrity each product platform is completely documented. Our internal configuration control will update changes immediately. Therefore all data and documentation to the component level are readily accessible.

Environmental Testing

ADLAS lasers are designed to withstand harsh environmental conditions. Each model is subjected to rigorous EMI, shock and vibration as well as temperature testing.

Production


All components and parts manufactured by outside vendors undergo critical incoming tests. Test procedures are fastidiously followed and results are recorded. Only parts meeting our requirements are stocked for use.


Specific assembly and handling routines help to protect sensitive components from damage through ESD, contamination, etc. during the manufacturing process. Critical components are documented on individual test records. Quality Circles and scheduled staff training courses are designed to raise the overall quality consciousness of our personnel. Assembly of critical components is conducted in our Clean Room facilities. Measurement instruments are under constant control and routine calibration.

Final Test

Each individual unit is tested to meet the key parameters of output power, polarization, beam quality, etc.. In addition, each laser is tested over an extended burn-in time under real operating conditions.

The persistently high level of quality awareness during development and manufacturing is the foundation of highly reliable, maintenance-free ADLAS lasers.

Milestones

- 1986: Incorporation of ADLAS Germany by Dr. Steyer and Dr. Kortz, both with many years of experience in development, production, and marketing of lasers.
- June 1986: world's first commercial diode pumped Nd:YAG laser product with 40 mW of TEM₀₀ output power.
- 1987: first commercial frequency doubled DPY with 2 mW of output power.
- 1988: first ADLAS q-switched laser
- 1989: optical output power IR/VIS increased to 350/80 mW.
- 1990: 750 mW/140 mW IR/VIS optical output power
- 1990: single longitudinal mode laser with 1 W output power and a linewidth of 10 kHz/10 ms.
- 1992: Incorporation of ADLAS, Inc. USA in Stow, MA
- 1993: 3000 mW/400 mW IR/VIS optical output power, more than 350 μJ q-switched pulse energy in the IR

Corporate Mission

ADLAS specializes in the development and production of diode laser pumped solid state lasers. ADLAS invests more than 20% of corporate revenues in the continued research and development of innovative pumping concepts and optimized laser configurations. ADLAS sees itself as a manufacturer of OEM laser light sources, capable of providing special lasers to fulfill the needs of our customers. ADLAS develops, assembles and tests all components in house. Active contact with both our customers and suppliers provides the flexibility to react to market requirements with optimally tailored products.

Corporate Goals

- Further miniaturization and power extension of our solid state laser products.
- Expansion of our leading market position as a supplier of diode pumped laser sources.
- Continually ensure the high level of our product quality.
- Adaption of our diode pumped solid state lasers to meet the requirements of new applications.

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■

430 Development Program


The design objectives of this program are to develop a compact laser with the following characteristics:

- 10mW CW at 430nm
- High wallplug efficiency
- Direct modulation capability
- Turn-key operation
- Zero maintenance package

To achieve this performance, we have been working with a single-frequency diode coupled to an external resonant doubler.

Coherent, Inc.

430 nm by Resonant Doubling of a Single Mode Diode

Coherent, Inc.

UPDATE—DPSS LASERS

New Advances in 430 nm, CW, Diode-Based Laser System

Coherent's Laser Group has been involved in the advancement of diode and diode-pumped laser technology for several years. These efforts have resulted in CW diode-pumped, Nd:YAG laser systems at 532 and 1064 nm. One of the current engineering programs is a CW, directly doubled diode (D^3) laser that produces 430 nm output.

Coherent's work has recently been supplemented by technology developed at IBM's Almaden Research Center. Combining the research efforts of both companies in the areas of diode stabilization and

externally resonant doubling has accelerated the commercial development of this technology.

D^3 Technology

To generate 430 nm, a single-frequency, 860 nm laser diode pumps an impedance-matched, external resonant cavity containing a KNbO_3 crystal. The external cavity must be stabilized to match the frequency of the input light to the cavity. This is accomplished using an FM locking technique. Greater than 10 mW CW output at 430 nm has been observed using a 100 mW, 860 nm laser diode.

System characterization and life testing are being conducted, and initial prototype systems are expected to be available later this year.


D^3 Applications

The D^3 has been a market-driven project since its inception. One of the major markets is information storage, whose

current systems are based largely on near infrared diodes. Switching to 430 nm will create a storage density 4 to ~6 times higher than currently possible.

Advantages also are anticipated for semiconductor inspection methods, where many commercial systems presently use air-cooled argon lasers operating at 488 nm. Because the resolving power of particle sizers is directly proportional to the laser wavelength, decreasing the wavelength of the light source will enhance the performance of these systems. An added benefit of D^3 technology is significantly reduced power consumption (20-50 times less than an air-cooled ion laser), which minimizes heat management requirements and lowers operation costs. A 430 nm laser is also of interest to the biotech community, for exciting fluorophors to study cellular metabolism.

For more information about these developments, please contact your local Coherent Sales Engineer or our Advanced Technical Sales Group at (800) 527-3786, or fax (800) 362-1170.

Optical diagram of D^3 .

Coherent, Inc., Laser Group, 5100 Patrick Henry Dr., Santa Clara, CA 95054 • Tel. 408-764-4983 • Fax 408-988-6838
E-mail: tech_sales@clg.com (Internet) • Benelux (079) 621313 • France (01) 6985 5145 • Germany (06074) 9140 •
United Kingdom (0223) 424065 • Japan (03) 3639-9811

266 Development Effort

OEM and scientific users have defined a need for a compact UV light source. We developed the following product definition for this source:

- 1-2mW CW at 266nm
- Wall plug efficient
- Turn-key operation
- Maintenance-free
- Thousands of hours of operation


The design that would best satisfy the product definition is a diode-pumped, frequency-quadrupled Nd:YAG laser.

Coherent, Inc.

266 Development Effort

The laser layout is shown below. We utilize an impedance-matched, external resonant doubler containing a BBO crystal pumped by a DPSS 532 laser.

We have observed >20mW CW output at 266nm with 120mW of 532nm pump power. The lifetime of this design is currently limited to ~100 hours. We are investigating methods to overcome this limitation.

Coherent, Inc.

12

UPDATE

Frequency-Quadrupled, Diode-Pumped, CW Nd:YAG Laser

UV Output Achieved with DPSS Series

Engineers at the Coherent Laser Group have demonstrated cw output at 266nm from a compact, solid-state laser. The results were obtained using an external enhancement cavity containing a BBO crystal pumped by a Coherent DPSS 532 laser. Conversion efficiencies of >20% were observed. The goal of this program is to develop a compact, reliable 1-2mW source at 266nm.

External Resonant Doubling

A single pass conversion efficiency of 0.02% $W^{-1} cm^{-1}$ is observed in BBO. Single pass conversion of 100mW of green in a 1cm crystal will

therefore generate ~2μW of UV. The use of an external resonant doubler can significantly improve the efficiency of UV generation. The design of the external cavity is shown in Figure 1. The cavity mirrors are high reflectors for the green output. The length of the external cavity is actively controlled to lock the cavity resonance to the frequency of the incident green light. This scheme results in a large enhancement of the green power and greatly increases the generated UV power, which scales as the square of the green intensity. Using resonant enhancement, several milliwatts of 266nm output was obtained by frequency doubling a Coherent DPSS 532 laser. Current experiments are investigating system stability and lifetime issues.

Applications

The development of a 266nm source has been application-driven. A com-

act UV laser will improve the sensitivity of forensic science techniques, because many natural and synthetic materials exhibit large UV absorbance. Researchers using confocal microscopes also have expressed a need for a UV source with a low M^2 value*, which would increase the system resolution compared to visible light sources. A 266nm laser will enhance the detection sensitivity of particle sizing devices, where the detection limit is directly proportional to the laser wavelength.

For more information about the 266 project or other DPSS lasers, contact Coherent's Advanced Technical Sales group at 1-800-527-3786 or fax 1-800-362-1170.

*For more information regarding the theory supporting the M^2 concept of measuring laser beam propagation characteristics, refer to "Beam Characteristics and Measurement of Propagation Attributes," by M.W. Sasnett and T.P. Johnston, Jr., SPIE Proceedings, vol. 1414, 1991.


Figure 1: Optical layout of the 266 prototype system.

Coherent Laser Group, 3210 Foster Drive, P.O. Box 10042, Palo Alto, CA 94303 • Telephone 415-855-7323 • Fax (408) 988-6838
France (01) 6985 5145 • Germany (06074) 9140 • United Kingdom (0223) 424065 • Japan (03) 3639-9811